


MICRO-ACCESS/MST INTRODUCER KITS AND SETS



GALT's Micro-Access/MST Introducing Kits and Sets were precisely designed for smaller access sites. Featuring 3cm, 5cm, 7cm, and 10cm lengths with the option of half-sized configurations. Each introducer features a smooth transition, easy and secure locking system, unique wing design, and funnel-shaped lead-in.

Features:

- Full kits and introducer sheath only
- 3cm, 5cm, 7cm, and 10cm lengths
- Half-sized configurations
- All kits include a Nitinol guidewire
- Safety components included



Smooth transition



Tearaway introducer available in half-sizes



For more information or to set up a meeting with a Galt Sales Consultant call (800) 639-2800.

MICRO-ACCESS/MST INTRODUCER KITS AND SETS

Micro-Access Tearaway Introducer Kits

Kit contents: PTFE Tearaway Introducer Set, a Nitinol Mandrel Guidewire, a safety guidewire introducer needle, a #11 safety scalpel, and a standard guidewire introducer needle

Kits contain a .018" x 40cm Marked Nitinol Mandrel Wire with a Stainless Steel Soft Tip

Catalog Number	Size	Description	Needle	Units/Box
KIT-010-00	3.0F	5cm	21ga x 5cm Safety	5
KIT-010-01	3.5F	5cm	21ga x 5cm Safety	5
KIT-010-02	4.0F	5cm	21ga x 5cm Safety	5
KIT-010-03	4.5F	5cm	21ga x 5cm Safety	5
KIT-010-04	5.0F	5cm	21ga x 5cm Safety	5
KIT-010-05	5.5F	5cm	21ga x 5cm Safety	5

Kits contain a .018" x 40cm Nitinol Mandrel Wire with a Stainless Steel Soft Tip

Catalog Number	Size	Description	Needle	Units/Box
KIT-010-37	6.5F	7cm	21ga x 7cm Safety	5

Kits contain a .018" x 65cm Marked Nitinol Mandrel Wire with a Stainless Steel Soft Tip

Catalog Number	Size	Description	Needle	Units/Box
KIT-010-12	4.0F	10cm	21ga x 7cm Safety	5
KIT-010-14	5.0F	10cm	21ga x 7cm Safety	5
KIT-010-15	5.5F	10cm	21ga x 7cm Safety	5
KIT-010-17	7.0F	10cm	21ga x 7cm Safety	5

Micro-Access Tearaway Introducer Sets

Set contents: PTFE Sheath and Dilator only

5cm Introducer Sets

Catalog Number	Introducer Size	Introducer Length	Accepts Wire Diameter	Units/Box
INT-106-00	3.0F	5cm	.018"	10
INT-106-01	3.5F	5cm	.018"	10
INT-106-02	4.0F	5cm	.018"	10
INT-106-03	4.5F	5cm	.018"	10
INT-106-04	5.0F	5cm	.018"	10
INT-106-05	5.5F	5cm	.018"	10
INT-106-06	6.0F	5cm	.018"	10

7cm Introducer Sets

Catalog Number	Introducer Size	Introducer Length	Accepts Wire Diameter	Units/Box
INT-107-30	3.0F	7cm	.018"	10
INT-107-40	4.0F	7cm	.018"	10
INT-107-50	5.0F	7cm	.018"	10
INT-107-60	6.0F	7cm	.018"	10

10cm Introducer Sets

Catalog Number	Introducer Size	Introducer Length	Accepts Wire Diameter	Units/Box
INT-106-10	3.0F	10cm	.018"	10
INT-106-12	4.0F	10cm	.018"	10
INT-106-13	4.5F	10cm	.018"	10
INT-106-14	5.0F	10cm	.018"	10
INT-106-15	5.5F	10cm	.018"	10
INT-106-16	6.0F	10cm	.018"	10
INT-106-17	7.0F	10cm	.018"	10

