

Biliary and Pseudocyst Drainage-Stents

Single Use



L-GES-SU-00-EN / LL-0017 Vers. 8.0 (01/2021)

Balanced Wall Diameters

Maximum inner lumen for form stability and kink resistance

Tapered and Round Endings

Atraumatic implantation

Pre-Assembled Stent Introducer-Sets

Convenient and time-saving positioning

Pseudo-Cyst Double-Pigtail Stents with Special Markers

Placement under EUS guidance



Medi-Globe GmbH


Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
Email: sales@medi-globe.de • www.medi-globe.de

Biliary Stent-Sets – Single Use

Biliary Stent-Sets – Single Use

Set Components based on Stent Diameter and Design:

	Stent-Diameter up to 7.0 French	Straight Stents with Diameter > 7.0 French	All Double-Pigtail Stents
Stent	✓	✓	✓
Guiding Catheter		✓	
Pusher	✓	✓	✓



Biliary Stent-Set – Straight Version

With bottom hole (Type „Amsterdam“); compatible with guide wire .035”

Catalogue No.	Stent		Guiding Catheter		Pusher	
	Size (French)	Length between Flaps (cm)	Size (French)	Length (cm)	Size (French)	Length (cm)
GBS-01-07-005S	7.0	5	n.a.	n.a.	7.0	220
GBS-01-07-007S	7.0	7	n.a.	n.a.	7.0	220
GBS-01-07-009S	7.0	9	n.a.	n.a.	7.0	220
GBS-01-07-012S	7.0	12	n.a.	n.a.	7.0	220
GBS-01-07-015S	7.0	15	n.a.	n.a.	7.0	220
GBS-01-08-005S	8.5	5	5.0	320	8.5	160
GBS-01-08-007S	8.5	7	5.0	320	8.5	160
GBS-01-08-009S	8.5	9	5.0	320	8.5	160
GBS-01-08-012S	8.5	12	5.0	320	8.5	160
GBS-01-08-015S	8.5	15	5.0	320	8.5	160
GBS-01-10-005S	10.0	5	6.0	320	10.0	160
GBS-01-10-007S	10.0	7	6.0	320	10.0	160
GBS-01-10-009S	10.0	9	6.0	320	10.0	160
GBS-01-10-012S	10.0	12	6.0	320	10.0	160
GBS-01-10-015S	10.0	15	6.0	320	10.0	160
GBS-01-11-005S	11.5	5	6.0	320	11.5	160
GBS-01-11-007S	11.5	7	6.0	320	11.5	160
GBS-01-11-009S	11.5	9	6.0	320	11.5	160
GBS-01-11-012S	11.5	12	6.0	320	11.5	160
GBS-01-11-015S	11.5	15	6.0	320	11.5	160




Biliary Stent-Set – Straight Version

With side ports; compatible with guide wire .035”

Catalogue No.	Stent		Guiding Catheter		Pusher	
	Size (French)	Length between Flaps (cm)	Size (French)	Length (cm)	Size (French)	Length (cm)
GBS-02-07-004S	7.0	4	n.a.	n.a.	7.0	220
GBS-02-07-006S	7.0	6	n.a.	n.a.	7.0	220
GBS-02-07-008S	7.0	8	n.a.	n.a.	7.0	220
GBS-02-07-011S	7.0	11	n.a.	n.a.	7.0	220
GBS-02-07-016S	7.0	16	n.a.	n.a.	7.0	220
GBS-02-08-004S	8.5	4	5.0	320	8.5	160
GBS-02-08-006S	8.5	6	5.0	320	8.5	160
GBS-02-08-008S	8.5	8	5.0	320	8.5	160
GBS-02-08-011S	8.5	11	5.0	320	8.5	160
GBS-02-08-016S	8.5	16	5.0	320	8.5	160
GBS-02-10-004S	10.0	4	6.0	320	10.0	160
GBS-02-10-006S	10.0	6	6.0	320	10.0	160
GBS-02-10-008S	10.0	8	6.0	320	10.0	160
GBS-02-10-011S	10.0	11	6.0	320	10.0	160
GBS-02-10-016S	10.0	16	6.0	320	10.0	160
GBS-02-11-004S	11.5	4	6.0	320	11.5	160
GBS-02-11-006S	11.5	6	6.0	320	11.5	160
GBS-02-11-008S	11.5	8	6.0	320	11.5	160
GBS-02-11-011S	11.5	11	6.0	320	11.5	160
GBS-02-11-016S	11.5	16	6.0	320	11.5	160

Biliary Stent-Sets – Single Use




Biliary Stent-Set – Double Pigtail

With position marker; compatible with guide wire .035"

Catalogue No.	Stent		Pusher	
	Size (French)	Length between Pigtails (cm)	Size (French)	Length (cm)
GBS-07-05-004S	5.0	4	5.0	220
GBS-07-05-006S	5.0	6	5.0	220
GBS-07-05-008S	5.0	8	5.0	220
GBS-07-05-009S	5.0	9	5.0	220
GBS-07-05-010S	5.0	10	5.0	220
GBS-07-05-012S	5.0	12	5.0	220
GBS-07-07-003S	7.0	3	7.0	220
GBS-07-07-004S	7.0	4	7.0	220
GBS-07-07-005S	7.0	5	7.0	220
GBS-07-07-006S	7.0	6	7.0	220
GBS-07-07-007S	7.0	7	7.0	220
GBS-07-07-008S	7.0	8	7.0	220
GBS-07-07-009S	7.0	9	7.0	220
GBS-07-07-010S	7.0	10	7.0	220
GBS-07-07-012S	7.0	12	7.0	220

Biliary Stent-Sets – Single Use




Biliary Stent-Set – Double Pigtail

With position marker; compatible with guide wire .035"

Catalogue No.	Stent		Pusher	
	Size (French)	Length between Pigtails (cm)	Size (French)	Length (cm)
GBS-07-08-003S	8.5	3	8.5	220
GBS-07-08-004S	8.5	4	8.5	220
GBS-07-08-005S	8.5	5	8.5	220
GBS-07-08-006S	8.5	6	8.5	220
GBS-07-08-007S	8.5	7	8.5	220
GBS-07-08-008S	8.5	8	8.5	220
GBS-07-08-009S	8.5	9	8.5	220
GBS-07-08-010S	8.5	10	8.5	220
GBS-07-08-012S	8.5	12	8.5	220
GBS-07-10-003S	10.0	3	10.0	220
GBS-07-10-004S	10.0	4	10.0	220
GBS-07-10-005S	10.0	5	10.0	220
GBS-07-10-006S	10.0	6	10.0	220
GBS-07-10-007S	10.0	7	10.0	220
GBS-07-10-008S	10.0	8	10.0	220
GBS-07-10-009S	10.0	9	10.0	220
GBS-07-10-010S	10.0	10	10.0	220
GBS-07-10-012S	10.0	12	10.0	220

Biliary Stent-Sets – Single Use


Biliary Stent-Sets – Single Use



Biliary Stent-Set – Double Pigtail for Common Bile Duct

Double pigtail with smaller diameter; with position marker


Catalogue No.	Stent		Pusher	
	Size (French)	Length between Pigtails (cm)	Size (French)	Length (cm)
GBS-17-07-004S	7.0	4	7.0	220
GBS-17-07-005S	7.0	5	7.0	220
GBS-17-07-006S	7.0	6	7.0	220
GBS-17-07-007S	7.0	7	7.0	220
GBS-17-07-008S	7.0	8	7.0	220
GBS-17-07-009S	7.0	9	7.0	220
GBS-17-07-010S	7.0	10	7.0	220
GBS-17-07-011S	7.0	11	7.0	220
GBS-17-07-012S	7.0	12	7.0	220



Pseudo-Cyst Stent-Set – Double Pigtail

Double Pigtail with Markers; compatible with guide wire .035"

Catalogue No.	Stent		Pusher	
	Size (French)	Length between Pigtails (cm)	Size (French)	Length (cm)
GBS-11-08-002S	8.5	2	8.5	220
GBS-11-08-005S	8.5	5	8.5	220
GBS-11-08-007S	8.5	7	8.5	220
GBS-11-08-009S	8.5	9	8.5	220
GBS-11-10-002S	10.0	2	10.0	220
GBS-11-10-005S	10.0	5	10.0	220
GBS-11-10-007S	10.0	7	10.0	220
GBS-11-10-009S	10.0	9	10.0	220




Biliary Stent-Set – Double Pigtail for Intrahepatic Bile Duct

Double pigtail with one smaller diameter; with position marker

Catalogue No.	Stent		Pusher	
	Size (French)	Length between Pigtails (cm)	Size (French)	Length (cm)
GBS-27-07-004S	7.0	4	7.0	220
GBS-27-07-005S	7.0	5	7.0	220
GBS-27-07-006S	7.0	6	7.0	220
GBS-27-07-007S	7.0	7	7.0	220
GBS-27-07-008S	7.0	8	7.0	220
GBS-27-07-009S	7.0	9	7.0	220
GBS-27-07-010S	7.0	10	7.0	220
GBS-27-07-011S	7.0	11	7.0	220
GBS-27-07-012S	7.0	12	7.0	220

Biliary Stents – Single Components – Single Use




Biliary Stent – Multi-Flap

Single components, not available in a set

Catalogue No.	Stent	
	Size (French)	Length between Flaps (cm)
GBS-04-07-005	7.0	5
GBS-04-07-007	7.0	7
GBS-04-07-009	7.0	9
GBS-04-07-012	7.0	12
GBS-04-07-015	7.0	15
GBS-04-08-005	8.5	5
GBS-04-08-007	8.5	7
GBS-04-08-009	8.5	9
GBS-04-08-012	8.5	12
GBS-04-08-015	8.5	15
GBS-04-10-005	10.0	5
GBS-04-10-007	10.0	7
GBS-04-10-008	10.0	8
GBS-04-10-009	10.0	9
GBS-04-10-010	10.0	10
GBS-04-10-012	10.0	12
GBS-04-10-015	10.0	15
GBS-04-11-005	11.5	5
GBS-04-11-007	11.5	7
GBS-04-11-009	11.5	9
GBS-04-11-012	11.5	12
GBS-04-11-015	11.5	15

Biliary Stents – Single Components – Single Use



Biliary Stent – Curved

Single components, not available in a set

Catalogue No.	Stent	
	Size (French)	Length between Flaps (cm)
GBS-05-07-005	7.0	5
GBS-05-07-007	7.0	7
GBS-05-07-009	7.0	9
GBS-05-07-012	7.0	12
GBS-05-07-015	7.0	15
GBS-05-08-005	8.5	5
GBS-05-08-007	8.5	7
GBS-05-08-009	8.5	9
GBS-05-08-012	8.5	12
GBS-05-08-015	8.5	15
GBS-05-10-005	10.0	5
GBS-05-10-007	10.0	7
GBS-05-10-009	10.0	9
GBS-05-10-012	10.0	12
GBS-05-10-015	10.0	15
GBS-05-11-005	11.5	5
GBS-05-11-007	11.5	7
GBS-05-11-009	11.5	9
GBS-05-11-012	11.5	12
GBS-05-11-015	11.5	15

Biliary Stents – Single Components – Single Use




Biliary Stent – Straight

Straight; with bottom hole (Type „Amsterdam“)

Catalogue No.	Stent	
	Size (French)	Length between Flaps (cm)
GBS-01-07-004	7.0	4
GBS-01-07-005	7.0	5
GBS-01-07-006	7.0	6
GBS-01-07-007	7.0	7
GBS-01-07-009	7.0	9
GBS-01-07-010	7.0	10
GBS-01-07-011	7.0	11
GBS-01-07-012	7.0	12
GBS-01-07-014	7.0	14
GBS-01-07-015	7.0	15
GBS-01-07-016	7.0	16
GBS-01-08-004	8.5	4
GBS-01-08-005	8.5	5
GBS-01-08-006	8.5	6
GBS-01-08-007	8.5	7
GBS-01-08-008	8.5	8
GBS-01-08-009	8.5	9
GBS-01-08-010	8.5	10
GBS-01-08-011	8.5	11
GBS-01-08-012	8.5	12
GBS-01-08-014	8.5	14
GBS-01-08-015	8.5	15

Catalogue No.	Stent	
	Size (French)	Length between Flaps (cm)
GBS-01-10-005	10.0	5
GBS-01-10-006	10.0	6
GBS-01-10-007	10.0	7
GBS-01-10-008	10.0	8
GBS-01-10-009	10.0	9
GBS-01-10-010	10.0	10
GBS-01-10-011	10.0	11
GBS-01-10-012	10.0	12
GBS-01-10-013	10.0	13
GBS-01-10-014	10.0	14
GBS-01-10-015	10.0	15
GBS-01-10-016	10.0	16
GBS-01-10-018	10.0	18
GBS-01-11-003	11.5	3
GBS-01-11-004	11.5	4
GBS-01-11-005	11.5	5
GBS-01-11-006	11.5	6
GBS-01-11-007	11.5	7
GBS-01-11-008	11.5	8
GBS-01-11-009	11.5	9
GBS-01-11-010	11.5	10
GBS-01-11-011	11.5	11
GBS-01-11-012	11.5	12
GBS-01-11-013	11.5	13
GBS-01-11-014	11.5	14
GBS-01-11-015	11.5	15

Biliary Stents – Single Components – Single Use




Biliary Stent – Straight

Straight; with side ports

Catalogue No.	Stent	
	Size (French)	Length between Flaps (cm)
GBS-02-05-004	5.0	4
GBS-02-05-005	5.0	5
GBS-02-07-004	7.0	4
GBS-02-07-006	7.0	6
GBS-02-07-008	7.0	8
GBS-02-07-010	7.0	10
GBS-02-07-011	7.0	11
GBS-02-07-013	7.0	13
GBS-02-07-016	7.0	16
GBS-02-08-004	8.5	4
GBS-02-08-006	8.5	6
GBS-02-08-008	8.5	8
GBS-02-08-011	8.5	11
GBS-02-08-013	8.5	13
GBS-02-08-016	8.5	16

Catalogue No.	Stent	
	Size (French)	Length between Flaps (cm)
GBS-02-10-004	10.0	4
GBS-02-10-006	10.0	6
GBS-02-10-008	10.0	8
GBS-02-10-009	10.0	9
GBS-02-10-010	10.0	10
GBS-02-10-011	10.0	11
GBS-02-10-012	10.0	12
GBS-02-10-013	10.0	13
GBS-02-10-016	10.0	16
GBS-02-11-004	11.5	4
GBS-02-11-006	11.5	6
GBS-02-11-008	11.5	8
GBS-02-11-011	11.5	11
GBS-02-11-013	11.5	13
GBS-02-11-016	11.5	16

Biliary Stents – Single Components – Single Use



Biliary Stent – Double Pigtail

Double pigtail with position marker


Catalogue No.	Stent	
	Size (French)	Length between Pigtails (cm)
GBS-07-05-004	5.0	4
GBS-07-05-006	5.0	6
GBS-07-05-008	5.0	8
GBS-07-05-009	5.0	9
GBS-07-05-010	5.0	10
GBS-07-05-012	5.0	12
GBS-07-07-003	7.0	3
GBS-07-07-004	7.0	4
GBS-07-07-005	7.0	5
GBS-07-07-006	7.0	6
GBS-07-07-007	7.0	7
GBS-07-07-008	7.0	8
GBS-07-07-009	7.0	9
GBS-07-07-010	7.0	10
GBS-07-07-012	7.0	12

Biliary Stent – Double Pigtail

Double pigtail with position marker

Catalogue No.	Stent	
	Size (French)	Length between Pigtails (cm)
GBS-07-08-002	8.5	2
GBS-07-08-003	8.5	3
GBS-07-08-004	8.5	4
GBS-07-08-005	8.5	5
GBS-07-08-006	8.5	6
GBS-07-08-007	8.5	7
GBS-07-08-008	8.5	8
GBS-07-08-009	8.5	9
GBS-07-08-010	8.5	10
GBS-07-08-011	8.5	11
GBS-07-08-012	8.5	12
GBS-07-08-014	8.5	14
GBS-07-08-015	8.5	15
GBS-07-10-002	10.0	2
GBS-07-10-003	10.0	3
GBS-07-10-004	10.0	4
GBS-07-10-005	10.0	5
GBS-07-10-006	10.0	6
GBS-07-10-007	10.0	7
GBS-07-10-008	10.0	8
GBS-07-10-009	10.0	9
GBS-07-10-010	10.0	10
GBS-07-10-012	10.0	12

Biliary Stents – Single Components – Single Use



Biliary Stent – Double Pigtail

Double pigtail with position marker

Catalogue No.	Stent	
	Size (French)	Length between Pigtails (cm)
GBS-17-07-004	7.0	4
GBS-17-07-005	7.0	5
GBS-17-07-006	7.0	6
GBS-17-07-007	7.0	7
GBS-17-07-008	7.0	8
GBS-17-07-009	7.0	9
GBS-17-07-010	7.0	10
GBS-17-07-011	7.0	11
GBS-17-07-012	7.0	12

Catalogue No.	Stent	
	Size (French)	Length between Pigtails (cm)
GBS-27-07-004	7.0	4
GBS-27-07-005	7.0	5
GBS-27-07-006	7.0	6
GBS-27-07-007	7.0	7
GBS-27-07-008	7.0	8
GBS-27-07-009	7.0	9
GBS-27-07-010	7.0	10
GBS-27-07-011	7.0	11
GBS-27-07-012	7.0	12



Pseudo-Cyst Stent – Double Pigtail

Pigtail with markers; compatible with guide wire .035"

Catalogue No.	Stent	
	Size (French)	Length between Pigtails (cm)
GBS-11-08-002	8.5	2
GBS-11-08-005	8.5	5
GBS-11-08-007	8.5	7
GBS-11-08-009	8.5	9
GBS-11-10-002	10.0	2
GBS-11-10-005	10.0	5
GBS-11-10-007	10.0	7
GBS-11-10-009	10.0	9

Stent Introducer Sets – Single Use

Stent Introducer Sets – Single Components – Single Use

Introducer Sets for Biliary Stents

- Combination of guiding catheter and pusher
- User-friendly and time-saving positioning of straight and curved biliary stents that require a guiding catheter
- With side injection port for contrast agent delivery over the horizontal guidewire and in every phase of stent placement
- Radiopaque metal markings support a controlled stent positioning under X-ray



Introducer Set for Biliary Stents

With 1 metal ring on the guiding catheter

Catalogue No.	Stent Size (French)	Pusher		Guiding Catheter	
		Length (cm)	Color Code	Size (French)	Length (cm)
GBS-31-08-180	8.5	180	Orange	5.0	220
GBS-31-10-180	10.0	180	Red	6.0	220
GBS-31-11-180	11.5	180	Blue	6.0	220



Introducer Set for Biliary Stents

With 4 metal rings at a distance of 5 cm on the guiding catheter

Catalogue No.	Stent Size (French)	Pusher		Guiding Catheter	
		Length (cm)	Color Code	Size (French)	Length (cm)
GBS-32-08-180	8.5	180	Orange	5.0	220
GBS-32-10-180	10.0	180	Red	6.0	220
GBS-32-11-180	11.5	180	Blue	6.0	220



Guiding Catheter

1 metal ring

Catalogue No.	Stent	
	Size (French)	Length (cm)
GBS-41-05-220	5.0	220
GBS-41-05-320	5.0	320
GBS-41-05-220	6.0	220
GBS-41-05-320	6.0	320



Guiding Catheter

4 metal rings at a distance of 5 cm

Catalogue No.	Stent	
	Size (French)	Length (cm)
GBS-43-06-220	6.0	220
GBS-43-06-320	6.0	320



Pusher

1 metal ring

Catalogue No.	Size (French)	Color Code	Length (cm)
GBS-42-05-220	5.0	Grey	220
GBS-42-07-160	7.0	Yellow	160
GBS-42-07-220	7.0	Yellow	220
GBS-42-08-160	8.5	Orange	160
GBS-42-08-180	8.5	Orange	180
GBS-42-08-220	8.5	Orange	220
GBS-42-10-160	10.0	Red	160
GBS-42-10-180	10.0	Red	180
GBS-42-10-220	10.0	Red	220
GBS-42-11-160	11.5	Blue	160
GBS-42-11-180	11.5	Blue	180
GBS-42-11-220	11.5	Blue	220